


## RoweSpiralLines / RoweTransfer-Sets



- pressure resistance up to 21 bar (305 psi)
- standard luer connectors
- visual drug flow

| Art.-No. | Material (tube) | Length (cm) | Tubing Ø (mm)              | Specifics   | Connection | PU ((pcs./box) / pal.) |
|----------|-----------------|-------------|----------------------------|---|------------|------------------------|
| A-7003   | PVC             | 150         | 1.65 x 3.15                | -   | fLL/mLL    | 50 / 3,600             |
| A-7009   | PVC             | 180         | 1.65 x 3.15                | -   | fLL/mLL    | 50 / 3,600             |
| A-7012   | PVC             | -           | 3.00 x 4.10<br>2.50 x 4.10 | RoweTransfer-Set<br>Tip: drip chamber<br>2 x Y-connector<br>2 x check valve | fLL/mLL    | 50 / 3,600             |
| A-7021   | PVC             | 180         | 1.65 x 3.15                | check valve,<br>Y-connector   | fLL/mLL    | 50 / 3,600             |
| A-7031   | PVC             | -           | 3.00 x 4.10<br>2.50 x 4.10 | RoweTransfer-Set<br>Tip: Spike<br>check valve                               | fLL/mLL    | 50 / 3,600             |
| A-7039   | PVC             | 150         | 1.65 x 3.15                | check valve,<br>Y-connector   | fLL/mLL    | 50 / 3,600             |
| A-7040   | PVC             | 210         | 1.65 x 3.15                | check valve,<br>Y-connector   | fLL/mLL    | 50 / 3,600             |
| A-7041   | PVC             | 210         | 1.65 x 3.15                | -   | fLL/mLL    | 50 / 3,600             |

customized systems are available

RoweMed AG - Medical 4 Life  
 www.rowemed.de • info@rowemed.de  
 Juri-Gagarin-Ring 4 • 19370 Parchim • Germany  
 Tel.: +49 (0) 3871 451280 • Fax : +49 (0) 3871 451282

