

RowePhil 18



- Hydrophilic injection filter for protection against particle and bacterial contamination of small amounts of drugs during administration
- Different connectors are available
- Filter area 2.3 cm²

- Filter diameter 18 mm
- Compatible with all Luer-Lock connectors
- Low priming volume

Art.-No.	Material (body)	Filling volume (ml)	Connection	Filter	PU ((pcs./box) / pal.)
A-6200	MABS	0.2	fLL/mL	0.2 µm nylon*	1,000 / 16,000
A-6202	MABS	0.2	fLL/mL	5.0 µm PET	1,000 / 16,000
A-6203	MABS	0.2	fLL/mL	0.2 µm PET	1,000 / 16,000
A-6204	MABS	0.2	fLL/mL	0.2 µm PES	1,000 / 16,000
A-6225	MABS	0.2	fLL/mLL	0.2 µm nylon*	1,000 / 16,000
A-6227	MABS	0.2	fLL/mLL	5.0 µm PET	1,000 / 16,000
A-6232	MABS	0.2	fLL/mLL	0.2 µm PES	1,000 / 16,000
A-6233	MABS	0.2	fLL/mLL	0.2 µm PET	1,000 / 16,000

* positively charged nylon membrane with endotoxin retention

RoweMed AG - Medical 4 Life
 www.rowemed.de • info@rowemed.de
 Juri-Gagarin-Ring 4 • 19370 Parchim • Germany
 Tel.: +49 (0) 3871 451280 • Fax : +49 (0) 3871 451282

