

## RoweFil



- Integrated venting through a PTFE membrane
- Compatible with all Luer-Lock connectors
- Can be connected to following systems easily
- Filter area 13.5 cm<sup>2</sup>
- Transparent connectors / tube / housing
- Visual drug flow control
- fLL/mLL
- MABS body

Art.-No.	Material (tube)	Length (cm)	Tubing Ø (mm)	Filling volume (ml)	Filter	Colour Code	OP (bar)	PU ((pcs./box) / pal.)
A-2356	PUR	10	1.5 x 2.7	2.20	0.2 µm nylon*	blue	3.2	50 / 5,600
A-2357	PUR	10	1.5 x 2.7	2.20	0.2 µm PET	light blue	3.2	50 / 5,600
A-2358	PUR	10	3.0 x 4.1	2.60	1.2 µm PET	yellow	2.5	50 / 5,600
A-2361	PUR	10	1.5 x 2.7	2.20	0.2 µm PES	-	3.2	50 / 5,600
A-2362	PUR	10 / 150	1.5 x 2.7	4.85	0.2 µm PES	-	3.2	240 / 3,840
A-2368	PUR	10	1.5 x 2.7	2.65	0.2 µm nylon * Y-part + clamp + rotating lock nut IN	blue	3.2	50 / 5,600
A-2373	PUR	2 x 10	3.0 x 4.1	3.40	1.2 µm PET	yellow	2.5	50 / 2,400
A-2374	PUR	10	1.5 x 2.7	2.40	0.2 µm nylon * check valve	blue	3.2	50 / 5,600
A-2377	PUR	10	1.5 x 2.7	2.65	0.2 µm nylon * Y-part + clamp + rotating lock nut HR	blue	3.2	50 / 5,600

\* positively charged nylon membrane with endotoxin retention for filtration up to 120 h

RoweMed AG - Medical 4 Life  
 www.rowemed.de • info@rowemed.de  
 Juri-Gagarin-Ring 4 • 19370 Parchim • Germany  
 Tel.: +49 (0) 3871 451280 • Fax : +49 (0) 3871 451282

